



Quantitative Analysis by XRD

Background

Full quantitative analysis is more difficult than the semi-quantitative work frequently performed with X-ray diffraction. Two methods for quantitative analysis are in common use - the Rietveld refinement method and the use of calibration standards. The Rietveld method is often termed the *Gold Standard* of quantitative analysis with accuracies in the 1% range, but use of calibration standards can give similar levels of accuracy. The difference, however, is that the Rietveld method is more widely applicable to a mixture of several phases, while the calibration standard method is limited to just 2 or 3 phases. Therefore, the Rietveld method is generally preferred as an all-purpose tool.

Example 1

Figure 1 shows a portion of the experimental diffraction pattern of a mixture of β -Tricalcium Phosphate (β -TCP) mixed with trace amounts of Hydroxyapatite (HA). Individual data points are shown as discrete dots and the computed profile is shown as the continuous line. At the top of the Figure is shown the *residual error* between the computed and experimental curves. Finally, at the bottom of the Figure are shown the stick patterns for the phases identified.

Applications

Corrosion products	Quantitative multiphase analysis
Forensic analysis	Amorphous/crystalline contents
Intermetallics	Quality control
Contaminants	Phase transformations
Mineral assays	Catalysts
Pharmaceuticals	Fiber analysis




Figure 1 Rietveld refinement of β -TCP with 1.7% HA

Once the phases have been identified and the necessary atomic information extracted from the database, a computed profile is calculated. Then, by use of a least-squares technique, the residual error is minimized by modifying the lattice parameters, the scale factors, atomic parameters such as the thermal vibrations, atomic positions, and occupancies. When the process is completed, the scale factors then provide the quantitative analysis.

Example 2

The method described in Example 1 is useful only when all of the details about the atomic and crystalline structure are known. In those cases where such information is not available, then the use of calibration standards is often used instead.

In this method, a number of control samples are made from the

pure materials over a composition range that overlaps the expected range in the target. Diffraction patterns are then prepared for each control sample, and the intensity ratio is prepared based upon two unobscured peaks from each phase. In the example shown in Figure 2, an X-ray diffraction pattern from the b-TCP/HA 1.5% mixture is shown.

One peak for each phase is highlighted in the Figure and a ratio of the intensities is computed. This process is then repeated for the other synthesized control samples. All of the data are used to construct the calibration curve shown in Figure 3. When an unknown sample is tested, the curve is used in reverse to estimate the quantitative amount through the measured intensity ratio.




Figure 2 XRD pattern of TCP/HA 1.5% mixture showing the non-overlapping peaks for the two phases used to generate the intensity ratio.




Figure 3 Calibration curve for TCP/HA mixtures

Our Pledge

H & M Analytical Services has over 30 years experience in X-ray diffraction. With our state-of-the-art equipment, we will strive to apply our experience and knowledge to solve your most challenging problems. In most cases, we will provide turnaround of 24 hours on phase identification analyses at no additional cost. Sample preparation services are also available. Most samples can be run on a flat fee basis. For details, please contact us.

Other Services We Provide

We provide a wide range of X-ray diffraction services, including:

- residual stress analysis
- precision lattice parameter
- texture analysis
- phase identification
- Rietveld analysis
- particle size determination
- high temperature XRD
- analysis of modulated films
- misfit strains
- fiber analysis
- crystal orientation
- grazing incidence angle
- retained austenite analysis



©H & M Analytical Services, Inc. 2003

H & M Analytical Services, Inc.
35 Hutchinson Road
Allentown, NJ 08501-1415
Tel: (609) 758-5700
Fax: (609) 758-5708
www.h-and-m-analytical.com